

ERWEKA DISINTEGRATION TESTER TYPE ZT222



On the ZT 222 series (available 1 to 4 test stations, operating simultaneously) the counting of the elapsed test duration is automatically started/ stopped when the arm with the basket is lowered/raised manually.

- ✓ 1 to 4 test racks
- ✓ 30 1 strokes/min.
- ✓ 55 2 mm stroke height
- ✓ Function keys and 2x LE display for entering nominal bath temperature and required run-time
- ✓ Starts/stops to count run-time automatically when basket rack is lowered
- ✓ Optional acoustic alarm on the completion of the pre-set run-time
- ✓ Manual lift-in and -out of the Baskets

DUTA MEDIKA

