

## MicroPette Plus Autoclavable Pipette Single-channel Adjustable Volume, Volume 2-10ml



Single-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume( $\mu$ l)	Maximуме Permissible Systematic error (Inaccuracy)		Maximуме Permissible Random error (Inaccuracy)	
			%	$\mu$ l	%	$\mu$ l
0.1-2.5 $\mu$ l	0.05 $\mu$ l	2.5	2.50%	0.0625	2.00%	0.05
		1.25	3.00%	0.0375	3.00%	0.0375
		0.25	12.00%	0.03	6.00%	0.015

8-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume ( $\mu$ l)	Maximуме Permissible Systematic error (Inaccuracy)		Maximуме Permissible Random error (Inaccuracy)	
			%	$\mu$ l	%	$\mu$ l
50-300 $\mu$ l	5 $\mu$ l	300	0.70%	2.1	0.25%	0.75
		150	1.00%	1.5	0.50%	0.75
		50	1.50%	0.75	0.80%	0.4

12-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume( $\mu$ l)	Maximуме Permissible Systematic error (Inaccuracy)		Maximуме Permissible Random error (Inaccuracy)	
			%	$\mu$ l	%	$\mu$ l
50-300 $\mu$ l	5 $\mu$ l	300	0.70%	2.1	0.25%	0.75
		150	1.00%	1.5	0.50%	0.75
		50	1.50%	0.75	0.80%	0.4

Fixed Volume Pipettes						
Volume Range	Increment	Test Volume( $\mu$ l)	Maximume Permissible Systematic error (Inaccuracy)		Maximume Permissible Random error (Inaccuracy)	
			%	$\mu$ l	%	$\mu$ l
5000 $\mu$ l	-	5000 $\mu$ l	0.3%	15	0.15%	7.5

