

Electrolyte Analyzer BG1000A



Features

- ✚ It applies advanced ISE (Ion-Selective Electrode) technology, convenient installation, easy maintenance, stable and reliable measurement.
- ✚ It has the advantages of automatic screensaver, friendly interface, and operate easily. Two working modes: Normal mode and saving models.
- ✚ It can save 1500 test results, and the data can be searched by date, ID, SN. The data will not be lost even if the power is cut off. Users can inquire or print at any moment.
- ✚ Revising data automatically to ensure the accuracy of data.
- ✚ It has a RS232 serial port, which can be connected to a computer. It may collect, manage and print measurement results by relevant software.

Specifications

Electrolyte Analyzer

- ✚ Sample: Serum, Whole Blood, Plasma, Diluted Urine
- ✚ Testing Time: not longer than 50s

Packaging Details: 355mm*350mm*420mm N.W.7kgs G.W.11kgs

